

**\*\*\*\*\* Medilux 2HF:**



Tube: 2x58W -direct Light.

Intensity: 3700lux

Size: 157x37x9cm

Weight: 13kg

**\*\*\*\*\* Medilux 2HFDO:**



Tube: 2x58W -direct Light.

Remote control: on/off

Intensity: 3700lux

Size: 157x37x9cm

Weight: 13kg

**\*\*\*\*\* Medilux 4HF:**



Tube: 4x58W -direct Light.

Intensity: 5000lux

Size: 157x67x9cm

Weight: 17kg

**\*\*\*\*\* Medilux 4HFDO:**



Tube: 4x58W -direct Light.

Remote control: on/off

Intensity: 5000lux

Size: 157x67x9cm

Weight: 17kg

**\*\*\*\*\* Medilux 4HFDO2:**



Tube: 4x58W -direct Light.

Remote control: on/off/2-intensity

Intensity: 5000lux/3000lux

Size: 157x67x9cm

Weight: 17kg

**\*\*\*\*\* Medilux 4HFDOST:**



Tube: 4x58W -direct Light.

Remote control: on/off/dimmable

Intensity: 5000lux~500lux

Size: 157x67x9cm

Weight: 17kg

**\*\*\*\*\* Medilux 4HF2:**



Tube: 4x58W -direct Light  
2x36W -indirect Light

Intensity: 5000lux  
Size: 157x67x9cm  
Weight: 18kg

\*\*\*\*\* **Medilux 4HF2DO:**



Tube: 4x58W -direct Light  
2x36W -indirect Light  
Remote control: on/off

Intensity: 5000lux  
Size: 157x67x9cm  
Weight: 18kg

\*\*\*\*\* **Medilux 6HF:**



Tube: 6x58W -direct Light.

Intensity: 5200lux  
Size: 157x97x9cm  
Weight: 24kg

\*\*\*\*\* **Medilux 6HFDO:**



Tube: 6x58W -direct Light.  
Remote control: on/off

Intensity: 5200lux  
Size: 157x97x9cm  
Weight: 24kg

\*\*\*\*\* **Medilux 6HFDO2:**



Tube: 6x58W -direct Light.  
Remote control: on/off/2-intensity

Intensity: 5200lux  
Size: 157x97x9cm  
Weight: 24kg


\*\*\*\*\* **Medilux 6HFDOST:**



Tube: 6x58W -direct Light.  
Remote control: on/off/dimmable

Intensity: 5200lux  
Size: 157x97x9cm  
Weight: 24kg

\*\*\*\*\* **Medilux 6HF2:**



Tube: 6x58W -direct Light  
2x36W -indirect Light  
Remote control: on/off/dimmable

**Intensity: 5200lux**  
Size: 157x97x9cm  
Weight: 24kg

**\*\*\*\*\* Medilux 6HF2DO:**



Tube: 6x58W -direct Light  
2x36W -indirect Light  
Remote control: on/off

**Intensity: 5200lux**  
Size: 157x97x9cm  
Weight: 24kg

**\*\*\*\*\* Remote control for Medilux:**



Dimmable version.

Battery: 12V